

Biological Indicator Tests

only \$1.79/test

BE COMPLIANT WITH CDC GUIDELINES



Your Complete Source for Sterility Assurance

1

Meet Mechanical Monitoring Guidelines

Per CDC:

"The mechanical monitors for steam sterilization include the daily assessment of cycle time and temperature by examining the temperature record chart (or computer printout) and an assessment of pressure via the pressure gauge."



2

Meet Chemical Monitoring Guidelines

Per CDC:

"Chemical indicators are affixed on the outside of each pack to show that the package has been processed through a sterilization cycle, but these indicators do not prove sterilization has been achieved. Preferably, a chemical indicator also should be placed on the inside of each pack to verify sterilant penetration. Chemical indicators usually are either heat-or chemical-sensitive inks that change color when one or more sterilization parameters (e.g., steam-time, temperature, and/or saturated steam; ETO-time, temperature, relative humidity and/or ETO concentration) are present."



3

Meet Biological Monitoring Guidelines

Per CDC:

"Steam and low temperature sterilizers (e.g., hydrogen peroxide gas plasma, peracetic acid) should be monitored at least weekly with the appropriate commercial preparation of spores. If a sterilizer is used frequently (e.g., several loads per day), daily use of biological indicators allows earlier discovery of equipment malfunctions or procedural errors and thus minimizes the extent of patient surveillance and product recall needed in the event of a positive biological indicator. Each load should be monitored if it contains implantable objects. If feasible, implantable items should not be used until the results of spore tests are known to be negative."



PREMIER MAIL-TO-LAB TESTING • ALL TYPES OF STERILIZERS • CALL FOR LOW PRICING
WEEKLY SPORE TESTING: \$164.00/Year • MONTHLY SPORE TESTING: \$74.00/Year

877-ATS-1199

www.SterilizerTesting.com

P.O. Box 270

Pearl River, NY 10965



Autoclave Testing Service, Inc.

Sterilizer Monitoring Since 1986

The Most Reliable, Convenient, and Affordable Method to Meet Guidelines and Regulations



UNIQUE TO ATS:
Includes plastic cover; helps to prevent dehydration of BIs. Digital read-out offers highly accurate temperature settings. Large capacity (24 ports) is great for busy offices!

RESULTS COMPLETED IN 24 HOURS!
TEST WEEKLY FOR ONLY \$1.79/TEST!

START-UP KIT INCLUDES:

Incubator (with 24 ports), ProTest BIs, 1 Box of 500 Class 5 Chemical Indicators, Log Record Book, BI Crusher

START-UP KITS

Start-Up Kit with 25 tests	\$295/kit
Start-Up Kit with 50 tests	\$335/kit
Start-Up Kit with 100 tests	\$375/kit

REFILLS

25 ProTest BIs	\$2.59 each
50 ProTest BIs	\$2.19 each
100 ProTest BIs	\$1.79 each

ProTest BI's:



The One Size Fits All Solution

**PROTEST IS COMPARABLE TO ALL OTHER BRANDS.
PROTEST IS 9MM WIDE AND INCUBATES BETWEEN 57-60°C.**

Call us for details or samples today!

ATS ALSO STOCKS OTHER STERILIZER MONITORING PRODUCTS:

Chemical Indicators • Indicator Tape • Sterilizer pouches & tubing • Bowie-Dick tests • BI Test Packs • Ultrasonic cleaners

877-ATS-1199

www.SterilizerTesting.com
P.O. Box 270
Pearl River, NY 10965



Autoclave Testing Service, Inc.
Sterilizer Monitoring Since 1986